

Software Development Kit

Get access to the SMI iView X™ eye tracking interfaces

- Integration of SMI iView X™ eye tracking systems with popular stimulus environments and custom applications
- API with high-level and low-level functions and sample code
- Access to real-time data and tracking status visualizations
- Built-in calibration and validation



Dr. Thomas Kieninger, German Research Center for Artificial Intelligence (DFKI):

“...the sample code of the SMI SDK was of great help in developing a custom application which integrates eye tracking data with object recognition technologies ...”

www.smivision.com

The SMI Software Development Kit

The SMI iView X™ Software Development Kit (SMI SDK) allows integration of the SMI iView X™ systems with popular stimulus software (e.g. MATLAB, PST E-Prime®, Python, NBS Presentation®) and custom applications written e.g. in C/C++ and .NET.

The SMI SDK enables these applications to easily remote-control the eye tracker and receive real-time data from the eye tracker by using the network-based communication interface of the iView X™ system.

High and low level functions with sample code

The SMI SDK contains an Application Programming Interface (API) consisting of several high-level functions (e.g. “calibration” and “tracking visualization”) and low-level functions (e.g. “get data”).

The extensive user manual provides a detailed description of more than 30 functions and ready-to-go sample code for a constantly increasing number of stimulus software products and programming languages.

Request the SMI SDK including extensive manual and sample code here: www.smivision.com/SDK

Example of C/C++ Sample Code

```
if (IV_Connect("127.0.0.1", 4444, "127.0.0.1", 5555) != 1)
    return;
if (IV_SetupCalibration(&calibrationSettings) != 1)
    return;
if (IV_Calibrate() != 1)
    return;
if (IV_GetSample(&sampleData) != 1)
    return;

cout << "Timestamp:" << sampleData.timestamp << endl;
cout << "Right Gaze X:" << sampleData.rightEye.gazeX << endl;
cout << "Right Gaze Y:" << sampleData.rightEye.gazeY << endl;
cout << "Right Diameter X:" << sampleData.rightEye.diam << endl;

IV_Disconnect();
```

Example of MATLAB Sample Code

```
loadlibrary ('View XAPI.dll', 'View XAPI.h');

pSample = ilbpointer ('SampleStruct', Sample);
pCalibrationSettings = ilbpointer ('CalibrationStruct', CalibrationSettings);

calllib ('View XAPI', 'IV_Connect', '127.0.0.1', int32 (4444), '127.0.0.1', int32 (5555))
calllib ('View XAPI', 'IV_SetupCalibration', pCalibrationSettings)
calllib ('View XAPI', 'IV_Calibrate')

calllib ('View XAPI', 'IV_GetSample', pSample);
get (pSample, 'Value')

calllib ('View XAPI', 'IV_Disconnect')

unloadlibrary ('View XAPI');
```

Remote-control and synchronization

The SMI SDK allows a stimulus software or custom application to remote-control and to synchronize with an iView X™ eye tracking system:

- Calibrate/validate the eye tracking system including visualization of the results
- Start/stop the experiment, and store the recorded eye data and separate trials
- Synchronize with the stimulus software through real-time messages

Access to real-time data and visualizations

The SMI SDK provides access to real-time eye data and online tracking status information, thus supporting the implementation of gaze-controlled applications and gaze-contingent paradigms:

- Real-time binocular gaze data, pupil diameter, online fixation events
- Low latency of eye data transfer
- Online tracking monitor and live view of the eye image (for selected devices)

Example of .NET Sample Code

```
if (IV_Connect("127.0.0.1".ToCharArray(), 4444, "127.0.0.1".ToCharArray(), 5555) != 1)
    return;

if (IV_SetupCalibration(ref calibrationSettings) != 1)
    return;
if (IV_Calibrate() != 1)
    return;
if (IV_GetSample(ref sampleData) != 1)
    return;

textBox1.Text = "Sample data - Timestamp:" + sampleData.Timestamp.ToString() +
    " - Right Gaze X:" + sampleData.rightEye.gazeX.ToString() +
    " - Right Gaze Y:" + sampleData.rightEye.gazeY.ToString() +
    " - Right Diameter X:" + sampleData.rightEye.diam.ToString() +
    IV_Disconnect();
```

Example of Python Sample Code

```
ViewXAPI = windll.LoadLibrary ("ViewXAPI.dll")

ViewXAPI.IV_Connect(c_char_p ("127.0.0.1"), c_int (4444), c_char_p ("127.0.0.1"), c_int (5555))

ViewXAPI.IV_SetupCalibration(byref(calibrationSettings))
ViewXAPI.IV_Calibrate()

ViewXAPI.IV_GetSample(byref (sampleData))

print "Timestamp: " + str(systemData.timestamp)
print "Right Gaze X: " + str(systemData.rightEye.gazeX)
print "Right Gaze Y: " + str(systemData.rightEye.gazeY)
print "Right Diameter X: " + str(systemData.rightEye.diam)

ViewXAPI.IV_Disconnect()
```

Technical Specification

Supported Operating Systems

- Windows XP 32/64 bit, Windows Vista 32/64 bit, Windows 7 32/64 bit, Linux/ Mac OS X (planned)

Sample Code

- C/C++, .NET, MATLAB, Python, PST E-Prime®, NBS Presentation®

Supported Eye Tracking Devices

- SMI RED, SMI RED250, SMI RED500, iView X™ HED, iView X™ Hi-Speed, iView X™ Hi-Speed Primate, iView X™ MRI, iView X™ TEM

For more information on SMI eye tracking solutions visit our website at www.smivision.com/egts.

Contact Information

SensoMotoric Instruments GmbH
Warthestr. 11
14513 Teltow
Germany
Phone: +49 (0) 3328 – 39 55 – 10
Fax: +49 (0) 3328 – 39 55 – 99
E-mail: sales@smi.de

SensoMotoric Instruments Inc.
28 Atlantic Ave
236 Lewis Wharf
Boston, MA 02110 USA
Phone: +1 - 617 - 557 - 00 10
Fax: +1 - 617 - 507 - 83 19
E-mail: sales@smivision.com



Request the SMI SDK!